KIRSCH BLOODHOUND XTN

XTREME TRACKING NOCK



Starter Kit includes Bloodhound XTN #1, two Practice Nocks, Tracker, Charger Kit, Spare Set of Bushings and Replacement Expandable Retention Mech, Magnetic Off Wand, Nocks.

Helping you recover more game with less grid searches, bring back your meat in less time with less spoilage, and less risk from other predators finding your trophy before you do.


- 1 Bushings (2) to snugly fit .244" (S), .246" (GT) arrow IDs, and .285", .297", and .300" bolt IDs.
- 2 Expandable Retention Mech. holds onto game.
- 3 Rechargeable battery w/ 1 year standby & 18-hours tracking.
- 4 Alum. front shoulder to accurately align to arrow/bolt end.
- 5 1.26" width Stopping Mech.
- 6 Embedded antennas.
- 7 Aircraft grade alum. core for high speed / high impact strength.
- 8 Symmetrical rear exposed LEDs (2) providing max viewing to shooter.
- 9 Charging pins (2) providing long reusable life.
- 10 Alum. rear shoulder for strength and accurate alignment with the nock.
- 11 Swappable pin nock for multiple crossbows (flat nock, halfmoon nock, and capture nock) and preferred arrow pin nock.



KIRSCH, LLC

P.O. Box 749, New Munster, WI 53152 800.970.8041

- + Lights super bright when shot, letting you know it is on.
- + Pulls off the back of your arrow during arrow pass through of the target.
- + Secures to the fatally shot game giving you a blood trail backup.
- + Interacts with your hand-held Tracker up to 400 yards away.
- + Tracker points you to any of up to 6 unique Bloodhound XTNs.
- Helps you recover your trophy when other methods fail.
- + Encouraging deadly arrow pass through.



www.KirschHunting.com CustomerService@KirschHunting.com

The hand held TRACKER

If you lose the bloodtrail, your backup is now in hand. Your personal Tracker uses RF (radio frequencies) communicating between the XTNs and Tracker for guiding you to your fatally shot game, recovering, not wasting, the land's bounty. No Bluetooth or cell service to worry about.



- 1 USB port for recharging battery
- 2 Integrated antenna doors (x4)
- 3 Power Button
- 4 Tracker battery life indicator
- 5 XTN battery life indicator
- 6 Directional arrow (to the XTN)
- 7 Bloodhound tracker wags tail when locked onto XTN
- 8 XTN # indicator, supports up to 6 unique XTNs
- 9 Distance indicator by %: 99% equals next to the XTN, 20% being possibly 400 yards away
- 10 Track Button to track a specific XTN #
- 11 Scan Button for finding any XTNs #s 1-6 within range
- 12 Stop Button to stop all action
- 13 Speaker Button (on/off)
- 14 Backlight Button (on/off)
- 15 Select Button
- 16 Loop bar for attaching wrist strap, etc.



KIRSCH, LLC P.O. Box 749, New Munster, WI 53152 800.970.8041

www.KirschHunting.com CustomerService@KirschHunting.com

Technical Specifications

Tracker

Housing Color: Black Weight: 10 ounces Dimensions (product): H: 6", W: 3.5", D: 1.5" Operating Temperature: -10'C to +50'C Storage Temperature: -20'C to +45'C up to 6m -20'C to +20'C up to 2yr Battery Type: Rechargeable Li-ion 24 hrs. in Search & Track Battery Life: Charge Cycles: 250 Range: Up to 400 yards Weatherproofing: Light Rain **XTN Tracking Nock** Wing Color: Translucent orange Metal Color: **Black Anodized** Translucent orange Pin Nock Color: LED Color: **Ultra-Bright Orange** Weight: XTN w/o pin nock: 95 grains Compound bow pin nock: 4 grains KIRSCH flat nock: 9.2 grains KIRSCH half-moon nock: 9.5 grains KIRSCH capture nock: 10.2 grains L: 3.4" total (w/o pin nock) Dimensions: L: 0.86" from rear of arrow shaft to front of pin nock W: 1.26" wing span Diameter: 0.345" max O.D. -30'C to +50'C Operating Temperature: Storage Temperature: -20'C to +45'C up to 6m -20'C to +20'C up to 2yr Battery Type: **Rechargeable Li-ion** Battery Life: 18 hrs. in Transmit mode 12 m in Standby & Storage Charge Cycles: 250 Weatherproofing: Water resistant, do not submerge Charging Temp: 0'C to +55'C Product Approvals: FCC part 15; IC

Limited Manufacturer's Warranty: Limited manufacturer's warranty guaranteed against defects in materials and workmanship not caused by normal wear, misuse, or improper maintenance as determined by Kirsch, LLC. It does not cover damage by impact from another arrow or impact with hard objects.



KIRSCH, LLC P.O. Box 749, New Munster, WI 53152 800.970.8041



Ordering Numbers – Compound Bow:

| Part Number | Use with Compound Bows shooting arrow Inside Diameters of .244" (S) & .246" (GT) |
|---|---|
| | Complete Kits (Bloodhound XTNs + Trackers) |
| | |
| XTN-31120 | Bloodhound XTN #1, 2 Practice XTNs, Tracker, wrist strap, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 3 nocks (Starter Kit) |
| XTN-31730 | Bloodhound XTNs #1-3, 3 Practice XTNs, Tracker, wrist strap, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks (Weekend Warrior Kit) |
| XTN-31960 | Bloodhound XTNs #1-6, 6 Practice XTNs, Tracker, wrist strap, 2 charger kits, magnetic off wand, spare set of bushings and replacement retention mechs, and 12 nocks (Professional Hunter Kit) |
| | Sets of XTNs (without Tracker) |
| | Bloodhound XTNs #2-3, 1 practice XTNs, magnetic off wand, spare set of bushings and |
| XTN-30A10 | replacement retention mechs, and 3 nocks (Expansion Kit A) |
| XTN-30730 | Bloodhound XTNs #1-3, 3 practice XTNs, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks |
| XTN-30830 | Bloodhound XTNs #4-6, 3 practice XTNs, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks (Expansion Kit B) |
| XTN-30960 | Bloodhound XTNs #1-6, 6 practice XTNs, 2 charger kits, magnetic off wand, spare set of bushings and replacement retention mechs, and 12 nocks |
| XTN-30030 | 3 practice XTNs, bushings, and 3 nocks |
| XTN-30060 | 6 practice XTNs, bushings, and 6 nocks |
| | Spares & Accessories |
| XTN-30100 | Bloodhound XTN #1, spare set of bushings and replacement retention mechs, magnetic off wand, and nock |
| XTN-30200 | Bloodhound XTN #2, spare set of bushings and replacement retention mechs, magnetic off wand, and nock |
| XTN-30300 | Bloodhound XTN #3, spare set of bushings and replacement retention mechs, magnetic off wand, and nock |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | Bloodhound XTN #4, spare set of bushings and replacement retention mechs, magnetic off wand, |
| XTN-30400 | and nock |
| XTN-30500 | Bloodhound XTN #5, spare set of bushings and replacement retention mechs, magnetic off wand, and nock |
| | Bloodhound XTN #6, spare set of bushings and replacement retention mechs, magnetic off wand, |
| XTN-30600 | and nock |
| XTN-30010 | Practice XTN - qty 1, spare bushing, and nock |
| XTN-30002 | Replacement nocks and bushings - 6 each |
| XTN-30003 | Replacement XTN springs and screws - 6 each |
| XTN-50001 | USB Charging Kit: USB charging adapter, USB cable, power supply |
| XTN-50004 | Bloodhound XTN magnet |



Ordering Numbers – Crossbow

| Part Number | Use with Crossbows shooting bolt Inside Diameters of .285", .297", & .300" and string inside angle at full draw of <78 degree w/ capture nock, <86 degree w/ flat nock & "V" nock |
|-------------|---|
| | Complete Kite (Pleadbound VTNa + Trackara) |
| | Complete Kits (Bloodhound XTNs + Trackers) Bloodhound XTN #1, 2 Practice XTNs, Tracker, wrist strap, charger kit, magnetic off wand, spare |
| XTN-41120 | set of bushings and replacement retention mechs, and 3 nocks of each: flat, "v", and capture (Starter Kit) |
| XTN-41730 | Bloodhound XTNs #1-3, 3 Practice XTNs, Tracker, wrist strap, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks of each: flat, "v", and capture(Weekend Warrior Kit) |
| XTN-41960 | Bloodhound XTNs #1-6, 6 Practice XTNs, Tracker, wrist strap, 2 charger kits, magnetic off wand, spare set of bushings and replacement retention mechs, and 12 nocks of each: flat, "v", and capture (Professional Hunter Kit) |
| | Sets of XTNs (without Tracker) |
| XTN-40A10 | Bloodhound XTNs #2-3, 1 practice XTNs, magnetic off wand, spare set of bushings and replacement retention mechs, and 3 nocks of each: flat, "v", and capture (Expansion Kit A) |
| XTN-40730 | Bloodhound XTNs #1-3, 3 practice XTNs, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks of each: flat, "v", and capture |
| XTN-40830 | Bloodhound XTNs #4-6, 3 practice XTNs, charger kit, magnetic off wand, spare set of bushings and replacement retention mechs, and 6 nocks of each: flat, "v", and capture (Expansion Kit B) |
| XTN-40960 | Bloodhound XTNs #1-6, 6 practice XTNs, 2 charger kits, magnetic off wand, spare set of bushings and replacement retention mechs, and 12 nocks of each: flat, "v", and capture |
| XTN-40030 | 3 practice XTNs, bushings, and 3 nocks of each: flat, "v", and capture |
| XTN-40060 | 6 practice XTNs, bushings, and 6 nocks of each: flat, "v", and capture |
| | Spares & Accessories |
| XTN-40100 | Bloodhound XTN #1, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40200 | Bloodhound XTN #2, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40300 | Bloodhound XTN #3, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40400 | Bloodhound XTN #4, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40500 | Bloodhound XTN #5, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40600 | Bloodhound XTN #6, spare set of bushings and replacement retention mechs, magnetic off wand, and nock of each: flat, "v", and capture |
| XTN-40010 | Practice XTN - qty 1, spare bushings, and nock of each: flat, "v", and capture |
| XTN-40002 | Replacement nocks and bushings - 6 of each: flat, "v", and capture |
| XTN-40003 | Replacement XTN springs and screws - 6 each |
| XTN-50001 | USB Charging Kit: USB charging adapter, USB cable, power supply |
| XTN-50004 | Bloodhound XTN magnet |

Errors, technical changes, and delivery delays may occur. Technical information refers only to the specifications of the products. Features may differ from the information provided.

